

NAME:


HONORS CHEMISTRY

SECTION:


LAB: ENDOTHERMIC AND EXOTHERMIC REACTIONS

Many chemical reactions give off energy. Chemical reactions that release energy are called *exothermic* reactions. Some chemical reactions absorb energy and are called *endothermic* reactions. You will study one exothermic and one endothermic reaction in this experiment.

In Part I, you will study the reaction between citric acid solution and baking soda. An equation for the reaction is:



In Part II, you will study the reaction between magnesium metal and hydrochloric acid. An equation for this reaction is:



Another objective of this experiment is for you to become familiar with using the DataMate data-collection program on the TI Graphing Calculator. In this experiment, you will use the program to collect and display data as a graph or list, to examine your experimental data values on a graph, and to print graphs and data lists.




Figure 1

Preparations

Problem

What types of energy changes occur in chemical reactions?

Objectives

- Observing energy changes in chemical reactions
- Using Vernier probes to collect and process data
- Using specific heat to calculate energy changes

Materials

LabPro or CBL 2 interface
TI Graphing Calculator
DataMate program
Temperature Probe
50-mL graduated cylinder
balance

Styrofoam cup
250-mL baker
citric acid, $\text{H}_3\text{C}_6\text{H}_5\text{O}_7$, solution
Sodium hydrogen carbonate, NaHCO_3
hydrochloric acid, HCl , solution
magnesium, Mg

Roles Record the names—make sure this is included in your lab report, too!

- Project Manager _____
Reads directions, keeps track of time, keeps group on task
- Quality Control Manager _____
Approves data table, checks data quality, checks calculations
- Materials Manager _____
Gets and returns materials, supervises sharing of materials

Safety

Wear goggles. Be careful not to use chipped or broken glassware. Hydrochloric acid solutions are caustic and may cause skin irritation. Wear safety goggles when working with these solutions. Immediately wash any spills or splashes on your skin. Call your teacher.

Do the Experiment

PROCEDURE

1. Copy the data table into your lab notebook. Be sure to record your observations and all relevant data.

Part I Citric Acid plus Baking Soda

2. Plug the Temperature Probe into Channel 1 of the LabPro or CBL 2 interface. Use the link cable to connect the TI Graphing Calculator to the interface. Firmly press in the cable ends.
3. Place a Styrofoam cup into a 250-mL beaker as shown in Figure 1. Measure out 30 mL of citric acid solution into the Styrofoam cup. Place the Temperature Probe into the citric acid solution.
4. Weigh out 10.0 g of solid sodium hydrogen carbonate on a piece of weighing paper.
5. Turn on the calculator and follow these steps to start the DATAMATE program.

TI-73 and TI-83 Calculators:

Press **PRGM**, then press the calculator key for the *number* that precedes the DATAMATE program (usually **1**). Press **ENTER**. You are now at the main screen of the program. Press **CLEAR** to reset the program.

TI-83 Plus Calculators:

Press **APPS**, then press the calculator key for the *number* that precedes the DATAMATE program. Press **ENTER**. You are now at the main screen of the program. Press **CLEAR** to reset the program.

TI-86 Calculators:

Press **PRGM**, press **F1** to select <NAMES>, and press a menu key to select <DATAM> (usually **F1**). Press **ENTER**. You are now at the main screen of the program. Press **CLEAR** to reset the program.

TI-89, TI-92, and TI-92 Plus Calculators:

Press **2nd** [VAR-LINK]. Use **▼** or the cursor pad to scroll down to "datamate", then press **ENTER**. Press **)** to complete the open parenthesis that follows "datamate" on the entry line and press **ENTER**. You are now at the main screen of the program. Press **CLEAR** to reset the program.

6. Set up the calculator and interface for the Temperature Probe.
 - a. Select SETUP from the main screen.
 - b. If the calculator displays a Temperature Probe in CH 1, proceed directly to Step 7. If it does not, continue with this step to set up your sensor manually.
 - c. Press **ENTER** to select CH 1.
 - d. Select TEMPERATURE from the SELECT SENSOR menu.
 - e. Select the Temperature Probe you are using (in °C) from the TEMPERATURE menu.
7. Set up the data-collection mode.
 - a. To select MODE, press **▲** once and press **ENTER**.
 - b. Select TIME GRAPH from the SELECT MODE menu.

- c. Select CHANGE TIME SETTINGS from the TIME GRAPH SETTINGS menu.
 - d. Enter "6" as the time between samples in seconds.
 - e. Enter "50" as the number of samples. The length of the data collection will be 5 minutes.
 - f. Select OK to return to the setup screen.
 - g. Select OK again to return to the main screen.
8. You are now ready to begin collecting data.
 - a. Select START on the main screen.
 - b. After about 20 seconds have elapsed, add the baking soda to the citric acid solution. Gently stir the solution with the Temperature Probe to ensure good mixing.
 - c. A real-time graph of temperature vs. time will be displayed on the calculator screen during data collection.
 - d. Temperature readings (in °C) can also be monitored in the upper-right corner of the graph.
 - e. Data collection will stop after 5 minutes, and a graph of temperature vs. time will be displayed.
 9. Dispose of the reaction products as directed by your instructor.
 10. Use the  or  keys (or the cursor pad on the TI-92) to examine the data points along the displayed curve of temperature vs. time. As you move the cursor right or left, the time (X) and temperature (Y) values of each data point are displayed below the graph. Determine the initial temperature, t_1 , and final (or minimum) temperature, t_2 . Record the temperature values in your data table (round to the nearest 0.1°C).
 11. Store the data from the first run so that it can be used later. To do this:
 - a. Press  to return to the main screen, then select TOOLS.
 - b. Select STORE LATEST RUN from the TOOLS MENU.

Part II Hydrochloric Acid Plus Magnesium

12. Measure out 30 mL of HCl solution into the Styrofoam cup. Place the Temperature Probe into the HCl solution. Note: The Temperature Probe must be in the HCl solution for at least 30 seconds before doing Step 15.
13. Obtain a piece of magnesium metal from your instructor.
14. Choose START on the main screen to begin data collection. After about 20 seconds have elapsed, add the Mg to the HCl solution. Gently stir the solution with the Temperature Probe to ensure good mixing. **Caution:** Do not breathe the vapors. Data collection will stop after 5 minutes.
15. Dispose of the reaction products as directed by your instructor. Rinse the Temperature Probe.
16. Examine the data points along the displayed curve of temperature vs. time. Determine the initial temperature, t_1 , and the final (or maximum) temperature, t_2 . Record the temperature values in your data table (round to the nearest 0.1°C).
17. You can also examine the data by viewing the data lists directly. To do this, press  to return to the main screen, and select QUIT. Then follow this procedure for your calculator:

TI-73, TI-83, and TI-83 Plus Calculators

To view the lists, press  to display the EDIT menu and then select Edit. Press  to scroll down through the data lists L1 (time) and L2 (temperature) for Part II. **Note:** When you choose to repeat a data collection, as you did in Step 11 of this experiment, the temperature data in L2 of the first data run will be stored in L3. You can view the data for Part I of the experiment in L1 and L3.

TI-86 Calculators

To view the lists, press $\boxed{2\text{nd}}$ [STAT] and select <EDIT> by pressing $\boxed{F2}$. Press $\boxed{\blacktriangledown}$ to scroll down through the data lists L1 (time) and L2 (temperature) for Part II. **Note:** When you choose to repeat a data collection, as you did in Step 11 of this experiment, the temperature data in L2 of the first data run will be stored in L3. You can view the data for Part I of the experiment in L1 and L3.

TI-89, TI-92, and TI-92 Plus Calculators

To view the data matrix, press $\boxed{\text{APPS}}$, select Data/Matrix Editor, then Current. Move the cursor down through the data lists c1 (time) and c2 (temperature) for Part II. **Note:** When you choose to repeat a data collection, as you did in Step 11 of this experiment, the temperature data of the first data run will be displayed in c3. You can view the data for Part I of the experiment in c1 and c3.

18. A good way to compare the curves is to view both sets of data on one graph.
 - a. Turn on the calculator and start the DATAMATE program.
 - b. Select GRAPH from the main screen, then press $\boxed{\text{ENTER}}$.
 - c. Select MORE, then select L2 AND L3 VS L1 from the MORE GRAPHS menu.
 - d. Both temperature runs should now be displayed on the same graph. Each point of Part I (citric acid and baking soda) is plotted with a box, and each point of Part II (hydrochloric acid and magnesium) is plotted with a dot.
19. (optional) Print a graph of temperature vs. time (with two curves displayed). Label each curve as "endothermic reaction" or "exothermic reaction".

DATA TABLE

	Part I	Part II
Final temperature, t_2	°C	°C
Initial temperature, t_1	°C	°C
Temperature change, Δt	°C	°C

Analyze and Apply

1. For each reaction, describe three ways you could tell a chemical reaction was taking place.
2. Which reaction took place at a greater rate? Explain your answer.
3. How much heat is lost as a 500. g cube of aluminum is cooled from 200.°C to 25.0 °C? The specific heat for aluminum is 0.897 J/g °C.
4. It takes 78.2 J to raise the temperature of 45.6 g lead by 13.3°C. Calculate the specific heat for lead.
5. A 28.2 g sample of nickel is heated to 99.8 °C and placed in a calorimeter containing 150.0 g water at 23.5 °C. After the metal cools, the final temperature of metal and water is 25.0 °C. Calculate the specific heat of nickel, assuming that all the heat lost by the nickel is gained by the water.

Address the following points in your discussion:

- Experimental evidence for exothermic vs. endothermic reactions
- Direction of energy flow relative to the system (the individual reactions) and the surroundings
- A summary of your findings
- Sources of error